

Sterile Adhesive Organizer

FAST ■ EASY ■ SECURE ■ VERSATILE



Self adhesive STRIP-T'S® hold and organize all surgical tubes, scopes and instruments any place on the surgical field.

- Economical - Easy to use
- Latex Free
- Saves pre-op time
- Versatile, can be used on every specialty surgery
- Several sizes and configurations
- Iridescent color for laser surgery
- Sterile shelf life 2 years

STRIP-T'S® is the least expensive apparatus to hold implements and is less expensive than tape, pouches or drapes. The STRIP-T'S® reduces the cost of replacing tubes, costly disposable suctions, surgical operative scopes and bovie, which tend to fall from the operative field.

CE 0050



**SURGICAL
CONCEPTS, INC.**


a division of KAPP SURGICAL INSTRUMENT, INC.

Worldwide Sales, Service & Manufacturing
4919 Warrensville Center Rd. Cleveland, Ohio 44128
Phone: 216-587-4400 / Fax: 216-587-0411 / 800-282-5277
www.kappsurgical.com / email: info@kappsurgical.com




FAST ■ EASY ■ SECURE ■ VERSATILE


Conveniently Holds and Organizes Surgical Implements In Every Surgical Procedure



LAPAROSCOPIC CHOLECYSTECTOMY



LAPAROSCOPIC CHOLECYSTECTOMY



GENERAL SURGERY
Cat. #SC-5000-3Y




ORTHOPEDIC SURGERY
Arthroscopy & total knee replacement shown with bovie scratch pad.



ORTHOPEDIC SURGERY
Total hip replacement.

STRIP-T'S® STANDARD CONFIGURATIONS



Cat. #SC-5000-3Y (Yellow) or SC-5000-3R (Red)

For use on all surgical procedures. Three self-closing Hook and Loop compartments plus one open Hook and Loop section. Two Hook and Loop pencil strips are provided to quickly wrap around your instrument or bovie and stick to the open Hook and Loop section. Scratch emery attached. Measures 2" x 12". Each Hook and Loop pencil strip measures 5/8" x 2-1/2".



Cat. #SC-5000-4Y (Yellow)

For general use on all surgeries. Four self-closing Hook and Loop compartments secure cords, tubes, scopes and instruments. Measures 2" x 10".



Cat. #SC-5000-4R (Red)

Note: Custom lengths and designs are available with 2,500 unit minimum order.



**SURGICAL
CONCEPTS, INC.**

a division of KAPP SURGICAL INSTRUMENT, INC.

Worldwide Sales, Service & Manufacturing
4919 Warrensville Center Rd. Cleveland, Ohio 44128
Phone: 216-587-4400 / Fax: 216-587-0411 / 800-282-5277
www.kappsurgical.com / email: info@kappsurgical.com